

Characterization of Thermal Radiative Properties of Space Materials at Cryogenic Temperatures

Tomáš Králík, Jiří Frolec

Institute of Scientific Instruments of the CAS, Brno, Czech Republic

Department of Magnetic Resonance and Cryogenics

Group of Cryogenics and Superconductivity

Workshop Academia & Industry in Space projects III 2026

Introduction

Thermal radiative properties for cryogenics and space applications

We can measure **total hemispherical emissivity or absorptivity 0.1 % - 100 %**

- Temperature of the radiator: **10 K - 320 K** (controlled)
- Temperature of the absorber: **5 K ~ 130 K** (sample depended)
- Sample dimensions: diameter **40 mm**
thickness max. **5 mm**
- Materials: solids, sheets, foils, coatings, thin films on one side of the sample


History

- First absorptivity measurements 2001-2002 for cryogenic applications
- Measurements for space applications since 2008

Collaboration

- Mostly contractual research
- Joint publication
- Consortium


JWST




MTG




Ariel



Athena



JUICE




Our experimental setup for thermal-radiative properties

We measure radiative heat exchange ($T_R \leftrightarrow T_A$) between two surfaces (\varnothing 40 mm disks)

Mutual emissivity $\epsilon_{RA} = \frac{Q_R}{A\sigma(T_R^4 - T_A^4)}$

Q_R ← Measured radiative heat flow [W]
 A ← area of the sample [m²]
 σ ← Stefan-Boltzmann constant [W·m⁻²·K⁻⁴]
 $T_R^4 - T_A^4$ ← Black body radiation



Apparatus in a Dewar vessel Scheme of the measurement chamber

Sample as radiator

Emissivity $\frac{1}{\epsilon_R} = \frac{1}{\epsilon_{RA}} - \frac{1}{\alpha_{REF}} + 1$

Sample as absorber


Absorptivity $\frac{1}{\alpha_A} = \frac{1}{\epsilon_{RA}} - \frac{1}{\epsilon_{REF}} + 1$

Sample in opposite position

Black reference sample $\epsilon_{REF}, \alpha_{REF} \approx 95\%$

Apparatus and samples


Opened measurement chamber



Black reference surface


Epoxy coating on Cu

Circular V grooved structure




Samples

Flat sample substrates



Sample substrates with screw



Solar white coating




Lunar regolith simulant




Examples of measurements

Emissivity and absorptivity of metals



Temperature of the radiator [K]

Emissivity of “black” coatings



Need more data? See our database or ask us.

Database of thermal radiative properties (ϵ , α) at cryogenic temperatures

Related paper: Frolec J., et al. 2019. Cryogenics. 97, 85-99

available at: <https://data.mendeley.com/datasets/z8t423rwwd/4>

- 69 measurements of about 50 individual samples
- temperature ranges 20 K-320 K; (plots & numeric data)
- mostly metallic materials with various surface finishes

